Gzunda GZS BED MOVER



The Gzunda GZS hospital bed mover allows one person to effortlessly move heavy hospital beds. The motorised unit operates equally well with or without a patient, simplifying the task of moving beds or patients between wards, departments, or to and from theatre, allowing the driver to focus on the needs of the patient during transport.

With a capacity of 600 kg, the silent motorised unit is able move a hospital bed up or down ramps, over carpet, into lifts, and around tight corners with minimal eff ort by the operator. The GZS has easy-to-use controls, including key start, twist grip variable speed control, horn, and emergency stop button. As with all Electrodrive equipment, the GZ10SL comes with a comprehensive 12-month warranty, including parts and labour.





Typical applications

Suitable for moving hospital beds. No under-bed clearance is required.

Features

Load capacity	Move up to 600 kg at 6° incline
Max. speed	Up to 5 km/hour
Speed mode	Three speed control with forward, reverse and emergency stop.
Usability	 Highly manoeuvrable in tight spaces. Intuitive controls with twist grip throttle and horn. No under-bed clearance required. No driver's licence required. Simple to use. Quiet, smooth operation. Zero emissions. Non-marking tyres.
Hitching	Push-button hitching—silent, fast and efficient.A variety of hitches are available.
Dimensions (L/W/H)	• GZS(IN): 1245/920/1230 mm • GZS(US): 1245/1015/1230 mm
Battery	Two 12V 26Ah MK-gel batteries with 24V smart charger.

Safety features

- Emergency stop button.
- Emergency back-off button.

ORDER CODE

Gzunda GZS 600 kg split tyne — no hitch (international)	GZSI(NH)
Gzunda GZS 600 kg split tyne — no hitch (US)	GZSUSI(NH)

